

# Motion Composer

Motion Composer is a software for collating, integrating and presenting interactive motion capture data. Embedded within Cortex, Motion Composer is able to continuously update information without exiting to a stand-alone software package, creating a seamless workflow. Users are able to create a user defined, graphical and interactive report comprising: Motion capture data, forceplate data, EMG data, kinetics and kinematics, documents, color video and live web content.

## Features

- Cortex plugin for presentation authoring
- Presentation tree/property controls for managing presentation resources
- Slide creation tools for laying out presentation slides
- Presentation loading/saving to presentation file format (XML)
- Presentation 'Pack-and-Go' for one-click packaging of viewer and data into a single file
- Freely redistributable viewer for presentation sharing New data views to support presentation views
- RTF editor pane for rtf/text create/import/view/edit
- HTML editor pane for HTML/media import/view/edit
- 2D graph pane for graphical information view/analysis

## Motion Composer Includes

### Motion Composer

Presentation authoring plug-in for Cortex. Author presentations seamlessly from your Cortex session.

### Motion View


Presentation player. Package and share presentations easily with this interactive viewer. Suitable for burning to a CD, emailing, or web distribution.

### New Data Views

New graph, document and multimedia views. Perfect for presentation data.

### Plug-in Kit

For advanced users, Visual Basic project templates for creating custom plug-in views.



Motion Analysis Corporation  
3617 Westwind Boulevard  
Santa Rosa, CA USA 95403

(T) 707.579.6500  
www.motionanalysis.com  
info@motionanalysis.com

